WITNESS® DIROFILARIA

Français	page	1
Español	página	5
Italiano	pagina	9
Português	página	13
English	page	17
Nederlands	bladzijde	21
Deutsch	Seite	25
Czech	strana	29
Slovak	strana	33
Hungarian	side	37
Bulgarian	страница	41

GÉNÉRALITÉS

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WITNESS® DIROFILARIA

GENERAL INFORMATION

Heartworm disease in dogs has a world wide distribution and is caused by the filarial nematode *Dirofilaria immitis*. Adult *D. immitis* worms are found in the right ventricle and pulmonary arteries of infected dogs. The adult female worm produces first stage larvae (microfilariae) which can circulate in the peripheral blood. When ingested by a suitable mosquito vector while taking a blood meal, these microfilariae undergo development to become infective third stage larvae. When these larvae are transferred to other dogs during feeding by these infected mosquitoes, they undergo further development for approximately 70 - 80 days.

After this time, immature adults begin to migrate to the distal pulmonary arteries and eventually grow to occupy the main pulmonary artery.

Some 6 - 8 months after infection by the third stage larvae, the parasites will have matured. Secreted adult heartworm antigen, as well as microfilariae, may then be found in the peripheral blood.

TEST INDICATION

An early and accurate heartworm diagnosis in an infected dog will allow treatment before the appearance of serious cardio-vascular troubles. Uninfected dogs can be protected from heartworm infection by using suitable prophylactic actions.

WITNESS® DIROFILARIA, test for the detection of the soluble antigen *D. immitis* is recommended for use when clinical history and/or clinical signs suggest heartworm disease or where a rapid verification of heartworm infection is required. The performance of the WITNESS® DIROFILARIA test is independent of circulating microfilariae and can be used to detect amicrofilaraemic or "occult" infections.

TEST PRINCIPLE

The WITNESS® DIROFILARIA test is a simple test, based on a Rapid Immuno Migration (RIM™) technology, using antibodies directed against specific epitopes of soluble adult *D. immitis* antigen.

Sensitized colloidal gold particles having bound to the antigen present within the sample (whole blood, serum or plasma) are allowed to migrate along a strip. The complex is then captured on a sensitized reaction zone where its accumulation causes the formation of a clearly visible pink/purple band. A control band, located at the end of the reading window, ensures that the test was performed correctly.

SPECIMEN INFORMATION

- The test can be performed on whole blood, serum or plasma. Whole blood samples must be anticoagulated with EDTA or heparin.
- Samples should always be collected with a sterile needle and syringe.
- Hemolysis does not significantly interfere with the test, but strongly hemolyzed samples may partly obscure (due to hemoglobin background) a weak positive result.

Storage

Sample should preferably be tested immediately after collection but not longer than 4 hours after collection, if stored at room temperature.

If testing is further delayed, samples should be refrigerated (+2°C - 8°C) for up to 2 days.

For prolonged storage, samples (serum and plasma only) should be kept frozen (-20°C).

KIT CONTENTS

- A. 10 pouches, each containing 1 test device and desiccant.
- B. 1 Buffer dropper bottle (2.8 ml).
- C. Instructions for use.
- D. 10 pipettes.

GENERAL PRECAUTIONS

1. Do not use components after expiration date.
2. Store the test kit at +2°C - 25°C. Do not freeze.
3. Use the test shortly after opening the sealed pouch (within 10 minutes).
4. Avoid touching or damaging membrane at windows (1), (2), (3).
5. The WITNESS[®] device should be placed on a flat, horizontal surface while performing the test.
6. Use a separate pipette for each sample.
7. Hold pipette and buffer bottle vertically when dispensing sample and buffer respectively.
8. Handle all reagents and samples as biohazardous material.
9. For veterinary use only.

TEST PROCEDURE AND RESULTS

Important : Allow samples and buffer drops to fall onto membrane at window (1). Do not touch pipette tip, sample or buffer drops, or buffer bottle tip directly to the membrane.

1. Sample application

- Tear open a pouch provided and place the test device on a flat horizontal surface for the duration of the test.
- Holding the provided pipette vertically, transfer one drop of sample to the sample well, window (1).

2. Buffer dispensing

- Check that the sample has truly absorbed into the membrane.
- Remove the cap from the buffer bottle, hold it vertically and add three drops of buffer to the sample well, window (1).
- Leave the test device flat during migration of sample/reagent complex through the reading window.

3. Reading test

- Wait 10 minutes, observe the presence or absence of pink/purple bands in reading windows (2) and (3).
- Sample results are read in window (2). The control band is read in window (3).

Note :

- *The test is complete and may be read before 10 minutes if two pink/purple bands are clearly visible in both windows (2) and (3).*
- *The presence of a pink/purple band only in window (3) before the 10 minutes does not mean that the test is complete. A pink/purple band in window (2) may develop slower than the control band in window (3).*

RESULTS

1. Validation

- Test is validated if a pink/purple band is present in the reading window (3).

2. Interpretation

- Negative for *D. immitis* antigen : No band in reading window (2), with one band in window (3).
- Positive for *D. immitis* antigen : One band in reading window (2), with one band in window (3).

Note :

- No band in control window (3) : invalid test.
- A test result should always be interpreted in the context of all available clinical information and history for the dog being tested.

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