# WITNESS® LEISHMANIA

<b>Français</b>	page	1
<b>Español</b>	página	5
<b>Italiano</b>	pagina	9
<b>Português</b>	página	13
<b>English</b>	page	17
<b>Nederlands</b>	bladzijde	21
<b>Deutsch</b>	Seite	25
<b>Czech</b>	strana	29
<b>Slovak</b>	strana	33

## GÉNÉRALITÉS

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# WITNESS® LEISHMANIA

## GENERAL INFORMATION

Phlebotomine sandflies vector or carry the disease agent associated with leishmaniasis. Sandflies can carry the disease agent from a variety of animals to humans. Although some forms of *Leishmania* are maintained by human populations, the majority of strains is more often associated with an animal reservoir and may not even infect humans (Bray 1974). Many forms of leishmaniasis are maintained by animal reservoirs in both the New and the Old World.

Three kinds of animal reservoirs can harbor *Leishmania* parasites (WHO Expert Committee 1990). A primary host harbors parasites and maintains the disease in the environment, acting as a stable reservoir. Because a primary host may have adapted to the parasite, the disease may not cause debilitation or illness in the host. The primary host plays the most important role in maintaining the epidemiology of the disease in animals (or humans) even if direct contact is not made. A secondary host may help maintain the parasite in the environment and bring the parasite into direct contact with animals (humans). A variety of wild and domestic animals can serve as *Leishmania* reservoirs. Rodents such as rats and gerbils as well as wild canines such as foxes are frequently implicated in outbreaks of the disease. Domestic dogs frequently harbor *L. infantum* which may be passed to humans, particularly children (Kirmse et al. 1987). Animals infected by parasites can show a variety of symptoms. Frequently, lesions develop on exposed skin (i.e., ears, nose, lips, tail) (Schlein et al. 1984). Less frequently, lesions develop on the neck or other parts of the body and result in severe depilation (Elbihart et al. 1987). Parasites can also cause problems in internal organs such as lungs, spleen, liver and kidneys (Duarte et al. 1986). Clinical signs of leishmaniasis in dogs can easily be confused with other diseases such as systemic lupus erythematosus, malignancies, and endocrine disorders (Harith et al. 1989). Some animals, however, show no obvious signs of the disease (Schlein et al. 1984 ; Who Expert Committee 1984 ; Jambou et al. 1986).

## TEST PRINCIPLE

The WITNESS® LEISHMANIA test is a simple test, based on Rapid Immuno Migration (RIM™) technology. Sensitized colloidal gold particles having bound to the antibody present within the sample (whole blood, serum or plasma) are allowed to migrate along a strip. The complex is then captured on a sensitized reaction zone where its accumulation causes the formation of a clearly visible pink band. A control band, located at the end of the reading window, ensures that the test was performed correctly.

## **SPECIMEN INFORMATION**

- The test can be performed on whole blood, serum or plasma. Whole blood samples must be anti-coagulated with citrate or heparin.
- Samples should always be collected with a sterile needle and syringe.
- Haemolysis does not significantly interfere with the test, but strongly haemolyzed samples may partly obscure (due to haemoglobin background) a weak positive result.

## **Storage**

Sample should preferably be tested immediately after collection but not longer than 4 hours after collection, if stored at room temperature.

If testing is further delayed, samples should be refrigerated (+2°C - 8°C) for up to 2 days.

For prolonged storage, samples (serum and plasma only) should be kept frozen (-20°C).

## **KIT CONTENTS**

- A. 10 pouches, each containing 1 test device and desiccant.
- B. 1 Buffer dropper bottle (2 ml).
- C. Instructions for use.
- D. 1 bag with 25 pipet tips.

### **Note :**

*Test and control bands are dyed yellow for quality control purposes. This dye does not interfere with the test results and will wash away while the test is developing.*

## **GENERAL PRECAUTIONS**

1. Do not use components after expiration date.
2. Store the test kit at +2°C - 25°C. Do not freeze.
3. Use the test shortly after opening the sealed pouch (within 10 minutes).
4. Avoid touching or damaging membrane at windows (1), (2), (3).
5. The WITNESS<sup>®</sup> device should be placed horizontally on a flat surface while performing the test.
6. Use a new pipette tip for each sample.
7. Hold pipette and buffer bottle vertically when dispensing sample and buffer respectively.
8. Handle all reagents and samples as biohazardous material.
9. For veterinary use only.

## TEST PROCEDURE AND RESULTS

*Important:* Allow samples and buffer drops to fall onto membrane at window (1). Avoid touching the membrane while applying the sample or buffer drops. Do not touch the membrane with pipette or buffer bottle tips.

### 1. Sample application

- Tear open a pouch provided and place the test device on a flat horizontal surface.
- Collect the sample with a pipette (not provided in the kit). Holding it vertically, transfer 2 times 5  $\mu\text{l}$  (or once 10  $\mu\text{l}$ ) of sample to the sample well (1).

### 2. Buffer dispensing

- Check that the sample has truly absorbed into the membrane.
- Remove the cap from the buffer bottle, hold it vertically and add four drops of buffer to the sample well (1).
- Leave the test device flat during migration of sample/reagent complex through the reading window.

### 3. Reading test

After 10 minutes, observe the presence or absence of pink bands in reading windows (2) and (3).

*Note :*

- *It is possible to read the test before 10 minutes if two pink bands (respectively in windows (2) and (3)) are clearly visible.*
- *The presence of only one band in reading window (3), prior to the end of the development time (10min), does not mean that the test is complete, as weak samples may appear slower than the control band.*

## RESULTS

### 1. Validation

Test is validated if a pink band is present in the reading window (3).

### 2. Interpretation

- No band in reading window (2), with one band in window (3): sample is negative for *Leishmania* antibodies.
- One band in reading window (2), with one band in window (3): sample is positive for *Leishmania* antibodies.

### Note :

- No band in control window (3): invalid test.
- A test result should always be interpreted in the context of all available clinical information and history for the dog being tested.
- It is recommended to re-test a whole blood sample as a serum or plasma sample, if the results do not correlate with the clinical symptoms.
- In case of doubtful results, sample should be re-tested after 6 or 8 weeks.

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