WITNESS® PARVO

Français	page	1
Español	página	5
Italiano	pagina	9
Português	página	13
English	page	17
Nederlands	bladzijde	21
Deutsch	Seite	25
Polish	strona	29
Hungarian	oldal	33
Swedish	sida	37
Danish	side	41
Finnish	sivu	45
Norwegian	siden	49
Czech	strana	53
Slovak	strana	57
Bulgarian	страница	61

INTRODUCTION

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INTRODUCTION

Canine Parvovirus (CPV) is member of the feline parvovirus subgroup. It is closely related to feline panleukopenia virus and mink enteritis virus, and is considered endemic to nearly all populations of domesticated and wild canines.

CPV causes two forms of disease : myocarditis and enteritis. Due to maternal antibody protection the myocardial form is rare. The enteric form, however, is prevalent and can be fatal to puppies and geriatric dogs. CPV enteritis causes severe, often bloody diarrhoea, vomiting, leukopenia and dehydration.

Transmission is faecal-oral and most infections occur from exposure to contaminated faeces. CPV is highly contagious and stable under a variety of environmental conditions. Rapid diagnosis of CPV allows for quarantine and prompt treatment of infected dogs. Diagnosis may be difficult in milder cases. WITNESS[®] PARVO detects all strains of Canine Parvovirus shed in faeces. Positive results with WITNESS[®] PARVO indicate the presence of canine parvovirus.

TEST PRINCIPLE

The WITNESS[®] PARVO test is a simple test, based Rapid Immuno Migration (RIM[™]) technology, which detects the presence of canine parvovirus antigen in canine faeces. Sensitized particles bound to the CPV antigen present within the sample (faecal extract) migrate along a membrane. The complex is then captured on a sensitized reaction zone where its accumulation causes the formation of a clearly visible purple band. A purple band in the control window ensures that the test was performed correctly.

SAMPLE INFORMATION AND STORAGE

Canine faecal material is required for this test. Stool samples may be stored at +2°C and 8°C for 48 hours. If longer storage is required the samples may be frozen at -20°C.

KIT CONTENTS

- A. 5 pouches, each containing 1 test device and desiccant.
- B. 5 disposable sample extraction pipettes.
- C. 1 extraction buffer bottle.
- D. Instructions for use.

PRECAUTIONS

- 1. Do not use components after expiration date.
- 2. Store the test kit at +2°C - 25°C. Do not freeze.
- 3. Use the test within 10 minutes after opening the sealed pouch.
- 4. Avoid touching or damaging the membrane at windows (1), (2) or (3).
- 5. The WITNESS® PARVO device should be placed on a flat horizontal surface while performing the test.
- 6. Use a separate extraction pipette for each sample.
- 7. Hold extraction pipette vertically when dispensing sample.
- 8. Handle all reagents and samples as biohazardous material.
- 9. For Veterinary Use Only.

"Healthy vaccinated animals may cause false positive results if tested from the 3rd to 10th day after vaccination."

TEST PROCEDURE AND RESULTS

1. SAMPLE COLLECTION

- Remove the top portion, containing the cotton swab, from a sample extraction pipette.
- Dispense Extraction Buffer into the base of the sample extraction pipette until the buffer reaches the graduation mark.
- Taking the top portion of the sample extraction pipette, coat the cotton swab with a thin layer of faeces from a stool sample, or an anal swab.
- Insert the swab into the dispensed extraction buffer in the base of the sample extraction pipette. Securely fit the top portion onto the base and tap to mix for 5 seconds to ensure good sample extraction.

2. SAMPLE DISPENSING

- Tear open the pouch provided and place the test device on a flat horizontal surface during the migration of the samples.
- Snap the top of the sample extraction pipette (at the blue line) containing the extracted faecal sample.
- Invert and hold the sample extraction pipette in a vertical (upright) position and squeeze the base gently to dispense five drops of sample to the sample well (Window (1)).
- If migration does not pass Window (2) by one minute, add one additional drop of sample in Window (1).

3. READING TEST

- Wait 10 minutes, then observe the presence or absence of purple bands in Windows (2) and (3).
- Sample results are read in Window (2). The control band is read in Window (3).

Notes :

- *The test is complete and may be read before 10 minutes if purple bands are visible in both Windows (2) and (3).*
- *The presence of a purple band only in Window (3) before 10 minutes does not mean that the test is complete. A purple band in Window (2) may develop slower than the control purple band in Window (3).*

4. RESULTS

Negative - A sample is negative for Canine Parvovirus antigen, if there is no purple band visible in Window (2) and a purple band in Window (3).

Positive - A sample is positive for Canine Parvovirus antigen if there is one purple band visible in Window (2) and one purple band in Window (3).

Note : A test is valid if a purple band is present in Window (3). If there is no purple band in Window (3), the test is invalid and should be repeated.

Remember : A test result should always be interpreted in the context of all available clinical information and history for the dog being tested.

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