# **WITNESS® FeLV**

<b>Français</b>	page	1
<b>English</b>	page	5
<b>Deutsch</b>	Seite	9
<b>Nederlands</b>	bladzijde	13
<b>Español</b>	página	17
<b>Italiano</b>	pagina	21
<b>Português</b>	página	25
<b>Polish</b>	strona	29
<b>Hungarian</b>	oldal	33
<b>Swedish</b>	sida	37
<b>Danish</b>	side	41
<b>Finnish</b>	sivu	45
<b>Norwegian</b>	siden	49
<b>Czech</b>	strana	53
<b>Slovak</b>	strana	57
<b>Bulgarian</b>	страница	61

## GENERALITES

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## GENERAL INFORMATION

Feline Leukaemia Virus (FeLV) is a contagious Retrovirus from the Mammalian Type-C group, which is present in many areas of the world. This presence is particularly endemic in high density or close contact populations.

Transmission occurs essentially by contact, mainly through saliva or nasal secretions and by biting or licking. Clinically healthy persistently viremic cats are known as a major source of the virus.

Vertical (congenital) transmission is also described.

Infection is characterised by the development of a viremia which can be followed by seroconversion with elimination of the pathogen. Infected cats may also present a chronic, persistent viremia which will lead to the development of both proliferative syndromes like lymphosarcoma or leukaemia and non-proliferative syndromes such as anaemia or immunosuppression, followed by death at short or middle term.

Diagnosis of FeLV infection is generally performed by the detection of viral antigen from the core protein, p27, which is produced in high quantities in viremic cats.

## TEST INDICATION

WITNESS® FeLV test is indicated for use when history and/or clinical signs may suggest infection by FeLV. It will also be particularly recommended prior to an FeLV vaccination, especially in cats belonging to an at risk population.

## TEST PRINCIPLE

The WITNESS® FeLV test is a simple test, based on Rapid Immuno Migration (RIM™) technology, using an antibody directed against an epitope of the p27 core protein.

Sensitised colloidal gold particles bound to p27 antigen present within the sample (whole blood, serum or plasma) are allowed to migrate along a strip. The complex is then captured on a sensitised reaction line where its accumulation causes the formation of a clearly visible pink band. A control band, located at the end of the reading window ensures that the test was performed correctly.

## SPECIMEN INFORMATION

- The test can be performed on unclotted whole blood anticoagulated with EDTA or heparin, serum or plasma.

- Samples should always be collected with a sterile needle and syringe.
- Haemolysis does not significantly interfere with the test, but strongly haemolysed samples may partly obscure a weak positive band.

### **Storage**

Anticoagulated whole blood samples should preferably be tested immediately after collection but not longer than 4 hours after collection, if stored at room temperature.

If testing is further delayed, samples should be kept refrigerated (+2°C and 8°C) for up to 4 days. For prolonged storage, samples (serum and plasma only) should be kept frozen (-20°C).

### **KIT CONTENTS**

- A. 10 pouches, each containing 1 test device and desiccant.
- B. 1 buffer dropper bottle (of 2.8 ml).
- C. Instructions for use.
- D. 10 pipettes.

### **GENERAL PRECAUTIONS**

1. Do not use components after expiration date.
2. Store the test kit at +2°C - 25°C. Do not freeze.
3. Use the test immediately after opening the sealed pouch. (within 10 minutes).
4. Avoid touching or damaging membrane at windows (1), (2), (3).
5. The WITNESS<sup>®</sup> device should be placed on a flat, horizontal surface while performing the test.
6. Use a separate pipette for each sample.
7. Hold pipette and buffer bottle vertically when dispensing sample and buffer respectively.
8. Handle all reagents and samples as biohazardous material.
9. For veterinary use only.

## TEST PROCEDURE AND RESULTS

*Important:* Allow samples and buffer drops to fall onto membrane at window (1). Do not touch pipette tip, sample or buffer drops, or buffer bottle tip directly to the membrane.

### 1. Sample application

- Tear open a pouch provided and place the test device on a flat horizontal surface.
- Holding the provided pipette vertically, transfer one drop of sample to the sample well (1).

### 2. Buffer dispensing

- Check that the sample has fully absorbed into the membrane.
- Remove the cap from the buffer bottle, hold it vertically and add three drops of buffer to the sample well (1).
- Leave the test device flat during migration of sample/reagent complex through the reading window.

### 3. Reading test

After 10 minutes, observe the presence or absence of pink bands in reading windows (2) and (3).

*Note:*

- *It is possible to read the test before 10 minutes if two pink bands (respectively in (2) and (3)) are clearly visible.*
- *The presence of only one band in reading window (3), prior to the end of the development time (10 mins), does not mean that the test is complete, as weak samples may appear slower than the control band.*

## RESULTS

### 1. Validation

Test is validated if a pink band is present in the reading window (3).

### 2. Interpretation

- No band in reading window (2), with one band in window (3): sample is negative for FeLV antigen.
- One band in reading window (2), with one band in window (3): sample is positive for FeLV antigen.

**Note:**

- No band in control window (3): invalid test.
- A test result should always be interpreted in the context of all available clinical information and history for the cat being tested.

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