

CRF



Canine Rheumatoid Factor Test Kit

Diagnose Canine Rheumatoid Arthritis Quickly and Accurately Before Severe Disease Progression

- **Rapid**
Results in 3 minutes
- **Simple**
1 Reagent and 2 Controls
- **Convenient**
No instrumentation needed
- **Flexible**
12-36 tests per kit

CRF

Canine Rheumatoid Factor Test Kit

CRF is the overwhelming choice of reference laboratories.

Simple. Rapid. Convenient.

For your laboratory, choose the diagnostic which provides you maximum value.

The Simple Procedure

1. Add sample & latex beads to slide.
2. Stir & then rotate for 3 minutes.
3. Read results.



The Disease

Canine Rheumatoid Arthritis (RA) is a potentially crippling disease that is difficult to diagnose. The disease is an erosive polyarthritis, frequently characterized by chronic and progressive multi-site lameness, joint swelling and joint destruction. Studies show that Rheumatoid Factor (RF) is present in 70% of dogs with rheumatoid arthritis.

The Diagnosis

Early diagnosis is complicated because, unlike human medicine, there is no patient input about symptoms.

Consequently, veterinary practitioners are faced with the problem of minimal or no information until the disease has progressed to a severe stage. By this time, the treatment for canine RA may be ineffective.

The Simple Solution

Synbiotics' simple and reliable CRF test kit provides a practical test methodology for canine RA factor. Using this test kit, one can quickly and routinely test for RF in dogs. Synbiotics' sensitive and specific immunoassay technique detects clinically significant levels of RF in less than five minutes.

Synbiotics developed the CRF test kit in response to the veterinarians need for rapid diagnostic information of this debilitating canine disorder. Now you can provide a test result within minutes which will support the veterinarians diagnosis.



12200 NW Ambassador Drive, Suite 101
Kansas City, Missouri USA 64163
1-800-228-4305
1-816-464-3500