


The 20-minute diagnostic for Equine Infectious Anemia



About the disease...

Equine Infectious Anemia is a viral disease affecting only horses and other Equidae. It is characterized by intermittent fever with viremia, anemia, progressive weakness and persistent infection for life. The virus apparently can be transmitted only by direct entry into the bloodstream with the help of biting insects, needles and other surgical instruments, or rarely through open wounds.

ViraCHEK[®]/EIA

- **Quick Answer:** The fastest EIA test available anywhere!
- **User-friendly:** Room temperature incubation and single wash step mean more convenience, less equipment, greater flexibility.
- **Easy to interpret:** Bright blue positives are read directly and visually.
- **Reliable:** Purified recombinant p26 antigen virtually eliminates non-specific reactions; provides greater confidence in test results.
- **18-month shelf life:** Extends usefulness of all kit components.
- **50µl sample size:** Conserves serum for other tests.
- **Manufactured by Synbiotics:** Your satisfaction is 100% guaranteed!

For further information or technical support for ViraCHEK[®]/EIA, contact Synbiotics Laboratory Sales at 800-841-1875.

©1998 Synbiotics Corporation SA/EIA-4/98-5M-L Printed in USA

SYNBIOTICS
CORPORATION

12200 NW Ambassador Drive, Suite 101
Kansas City MO 64163
Tel: ++816-464-3500, Fax: ++816-464-3521
www.synbiotics.com