



SYNBIOTICS
CORPORATION

ANTIGENE SAW

**CONCENTRATED ANTIGEN
FOR THE SEROLOGICAL DIAGNOSIS
OF BRUCELLOSIS BY SLOW AGGLUTINATION
(WRIGHT'S SERO-AGGLUTINATION)**

COMPOSITION

Concentrated suspension of *Brucella abortus* (Weybridge strain 99) inactivated by heat and phenol (0.5 %) and dispersed in a 0.5 % phenol buffer.

INDICATIONS

Serological diagnosis of bovine brucellosis by slow agglutination.

ADMINISTRATION

Shake well before use.

Priorly dilute to 1/10th in a 0.5 % phenol physiological solution.

STORAGE

Store between +2°C and +8°C, protected from light. Do not freeze.

PRESENTATION

100-ml bottle

CODE ASAW

ANTIGENE SAW
400 reactions

LEAFLET REVISED : July 2006

FOR IN VITRO USE

SYNBIOTICS EUROPE - 2 rue Alexander Fleming 69367 LYON Cedex 07 - FRANCE