



# FRESH EXPRESS® KIT

## *Inventory List*

*The following is a list of the items contained in the FRESH EXPRESS® Kit. Please check all the contents upon receipt of the kit. If any materials are missing, please contact Synbiotics Technical Service immediately at (800) 228-4305 or (816) 464-3500.*

- 1 FRESH EXPRESS® shipping box with Styrofoam® liner
- 2 Chill packs – ***Place chill packs in freezer upon arrival, for at least 24 hours.***
- 1 Vial of FRESH EXPRESS® Semen Extender – ***Refer to instructions for storage details.***
- 1 Vial of Semen Separating Solution – ***Refer to instructions for storage details.***
- 1 Small Styrofoam® box with cardboard sleeve and plastic zipper bag
- 1 FRESH EXPRESS® Canine Semen Shipping Tube
- 1 FRESH EXPRESS® Canine Insemination Syringe
- 1 Eight-inch insemination pipette
- 1 Twelve-inch insemination pipette
- 1 Semen Collection Evaluation and Insemination Form
- 1 FRESH EXPRESS® Extender Instructions
- 1 Semen Separating Solution Instructions
- 1 Collection/Insemination Record
- 1 Instructions for the Insemination
- 1 AKC Application to Register a Litter Resulting from Artificial Insemination Using Fresh Extended Semen

**Semen collection supplies are not contained in this kit and must be ordered separately. Newspaper is also recommended to wrap around the inner Styrofoam® box for added insulation.**

### **THESE ITEMS ARE NON-RETURNABLE**

**Note:** FRESH EXPRESS® Semen Extender is packaged in 10 ml aseptic vials. Extender can be stored in the refrigerator for up to 30 days. If longer storage is required, it can be stored in the freezer until the expiration date stated on the bottle (do not thaw and refreeze). Please see FRESH EXPRESS® Extender Instructions for more information.

Semen Separating Solution should **NOT** be frozen.